

Protocol MSU_MSMC_005 version 1.1

Preparation of internal standard cocktail solutions for analysis of animal tissue metabolites

Last updated: May 6, 2019 by A. D. Jones

Reagents

Compound name	Source	MSMC Standard Bar Code
Palmitic acid-d ₃₁	Sigma-Aldrich #366897	IS_0003
[¹³ C ₆]D-glucose	Sigma-Aldrich #389374	
L-Phenylalanine-d ₇	Toronto Research Chemicals #P319416	IS_0005
Pyruvic acid- ¹³ C ₃ sodium salt	Toronto Research Chemicals #P998904	IS_0006
(-)-Inosine- 1',2',3',4',5'- ¹³ C ₅	Toronto Research Chemicals #I661003	IS_0004

Materials

15-mL polypropylene centrifuge tubes

2-mL polypropylene microcentrifuge tubes

Pyrex screw cap culture tubes with PTFE-lined phenolic caps, 13 x 100 mm (Fisher Scientific #14-933A)

Glass culture tubes with PTFE-lined screw caps, 10 x 75 mm (BMB Stores #1496215)

Glass media bottle (25 or 125 mL)

Graduated cylinders

Solvents

HPLC grade acetonitrile

MilliQ water (> 18 MΩ)

HPLC grade chloroform

HPLC grade methanol

Prepare at least 100 mL of a solution of 50% acetonitrile (v/v) in MilliQ water.

Prepare at least 10 mL of a solution of 2:1 chloroform/methanol in a glass bottle with PTFE-lined cap

Prepare stock solutions as follows:

Stock 1	Dissolve 5 mg of palmitic acid-d ₃₁ in 5000 μL (preferably measured using a positive-displacement pipetter) of 2:1 chloroform:methanol in a 13 x 100 mm glass culture tube (that has a PTFE-lined cap). Cap the tube and vortex for 30 seconds. Label this as Stock #1 (1.0 mg/mL palmitic acid-d ₃₁). Wait at least 15
----------------	--

	minutes before withdrawing aliquots. Store at -20°C. Allow to thaw to room temperature for 1 hour before withdrawing liquid.
Stock 2	Dissolve 4 mg of [¹³ C ₆]D-glucose in 1000 µL of 50% acetonitrile in water in a 2-mL polypropylene microcentrifuge tube; Vortex until all of the solid has dissolved. Label this as Stock #2 (4 mg/mL [¹³ C ₆]D-glucose)
Stock 3	Dissolve 5 mg of phenylalanine-d ₇ in 1000 µL of 50% acetonitrile in water in a 2-mL polypropylene microcentrifuge tube; Vortex until all of the solid has dissolved. Label this as Stock #3 (5 mg/mL phenylalanine-d ₇).
Stock 4	Dissolve 5 mg of [¹³ C ₃]pyruvic acid sodium salt in 1000 µL of 50% acetonitrile in water in a 2-mL polypropylene microcentrifuge tube; Vortex until all of the solid has dissolved. Label this as Stock #4 (5 mg/mL [¹³ C ₃]pyruvic acid sodium salt)
Stock 5	Dissolve 1 mg of inosine-1',2',3',4',5'- ¹³ C ₅ in 1000 µL of 50% acetonitrile in water in a 2-mL polypropylene microcentrifuge tube; Vortex until all of the solid has dissolved. Label this as Stock #5 (1 mg/mL inosine-1',2',3',4',5'- ¹³ C ₅)

Prepare a mixed stock solution of the polar water-soluble internal standards by mixing 100. µL each of Stock 2, Stock 3, Stock 4, and Stock 5 in a 2-mL polypropylene microcentrifuge tube. Add an additional 600 µL of 50% acetonitrile/50% water, and vortex to mix. Label the tube as "Mixed polar internal standard cocktail". Concentrations will be:

Polar compounds

[¹³ C ₆]D-glucose	400 µg/mL	(10 µg in 25 µL)
Phenylalanine-d ₇	500 µg/mL	(12.5 µg in 25 µL)
[¹³ C ₃]pyruvic acid sodium salt	500 µg/mL	(12.5 µg in 25 µL)
inosine-1',2',3',4',5'- ¹³ C ₅	100 µg/mL	(2.5 µg in 25 µL)

Mix 500 µL of the palmitic acid-d₃₁ stock solution with 500 µL of 2:1 chloroform/methanol in an amber glass autosampler vial with PTFE-lined cap.

Nonpolar compounds

Palmitic acid-d ₃₁	500 µg/mL	(12.5 µg in 25 µL)
-------------------------------	-----------	--------------------