

## Protocol MSU\_MSMC\_011

### Extraction of lipids and conversion to fatty acid methyl esters (FAMES)

Purpose: Simultaneous extraction and transesterification of total fatty acids and detection by GC/MS.

Last modified: July 7<sup>th</sup>, 2019

#### Reagents:

Methanol - Sigma, Cat# 646377

Chloroform – VWR, Cat# MK4443 (06 for 1 L, 10 for 4 L bottle)

Butylated hydroxytoluene (BHT) – Sigma, Cat# W218405 (Also available in the MSMC)

Formic acid, 98% – Sigma, Cat#FX0440

Potassium chloride, KCl

Sodium chloride, NaCl

Pentadecanoic acid (C15:0), analytical standard – Sigma, Cat# 91446; an alternative (for milk, animal, or microbial samples) is [<sup>13</sup>C16]palmitic acid (MSMC standard IS\_0015)

3M Methanolic hydrogen chloride - Sigma, Cat# 90964

n-Hexane, Spartan marketplace, Cat#15541750 (hexane can also be replaced with isooctane)

#### Materials:

2.0 mL Eppendorf Safe-Lock tubes

3 mm Chrome Stainless steel ball bearings (#Kit12064, VXB.com) or suitable substitute (may be available from Mass Spec Core)

Borosilicate glass (pyrex) tubes: 13x100 mm or another tube of sufficient volume - Fisher Scientific, Cat# 14-933A

Phenolic caps with PTFE liners for glass tubes – Sigma, Cat# CLS999813

Glass autosampler vials, caps and limited volume inserts (available in the BMB store)

#### Equipment:

Pipettor (1000 µL and 200 µL) and pipet tips

Lab shaker or vortexer

Tabletop centrifuge

Speed vac or nitrogen gas evaporator

Water bath

Centrifuge with a rotor that will fit the pyrex glass tubes

#### Notes:

Recommendations for tissue amount to be extracted:

- Plant tissue (leaf): 10-50 mg of fresh tissue.
- Animal tissue: 10-50 mg
- Bacterial cells: 10-50 mg of fresh cells

#### Folch, Lee and Stanley extraction protocol (chloroform/methanol):

1. Prepare the extraction solvent of 2:1 chloroform/methanol containing 0.01% (w/v) BHT, 1% (v/v) formic acid and pentadecanoic acid (10 µg/mL) as an internal standard. Alternatively, [<sup>13</sup>C16] palmitic acid (also 10 µg/mL) may be added in place of pentadecanoic acid for samples where pentadecanoic acid may be expected as a constituent of the lipids.

2. Prepare a 0.88% (w/v) solution of potassium chloride in water (8.8 g/L).
3. Add one 3-mm ball bearing per Eppendorf tube and homogenize frozen sample using a lab shaker or vortexer. Spin down briefly to collect the sample at the bottom of the tube. Make sure tissue stays frozen during the homogenization process.
4. Add 1000  $\mu\text{L}$  of extraction solvent (from Step 1) and cap the tube. This will add 10  $\mu\text{g}$  of internal standard to the tube. Vortex thoroughly for 3 mins.
5. Add 333  $\mu\text{L}$  of 0.88% potassium chloride in water and cap the tube. Vortex carefully and thoroughly.
6. Centrifuge for 10 min at 13500  $\times g$  in a tabletop centrifuge at ambient temperature. There should be two phases present after centrifugation. The upper layer is the water/methanol fraction while the lower will be the chloroform/methanol fraction. The lipids should be present in the lower layer.
7. Carefully remove a fixed volume (e.g. 500  $\mu\text{L}$ ) of the lower chloroform/methanol layer to a pyrex glass tube. Be careful in not removing any of the upper aqueous layer.
8. Dry the chloroform/methanol to dryness using either a speed-vac or nitrogen gas evaporator (both are available in the MSMC facility)
9. The dried lipid extracts can be stored sealed at  $-80\text{ }^{\circ}\text{C}$  until the trans-esterification reaction is ready to be performed.

#### **Trans-esterification protocol:**

1. Prepare 1M methylation reagent. For example, mix 20 mL methanol with 10 mL of 3 M methanolic hydrogen chloride. You will need 1 mL for each sample to be transesterified.
2. Prepare 0.9% (w/v) sodium chloride in water (9 g/L).
3. Add 1000  $\mu\text{L}$  of the methylation reagent (1 M methanolic HCl) to each dried lipid extract (from above). Ensure that the tube is capped and sealed.
4. Heat tubes in an  $80^{\circ}\text{C}$  water bath for 1 hour (note, this is above the boiling point of methanol). Check the tubes after 5 min to ensure that the caps are sealing well and that solvent is not evaporating during the trans-esterification reaction. It is often helpful to open the tube briefly to allow air to escape before re-sealing the cap.
5. Remove the tubes from the water bath. Let them cool to room temperature before opening.
6. Add 150  $\mu\text{L}$  hexane (or isooctane) and 1000  $\mu\text{L}$  0.9% (w/v) sodium chloride in water. Vortex thoroughly for 1 min.
7. Spin 10 min at 1,500  $\times g$  in a centrifuge with an appropriate rotor for pyrex glass tubes.
8. Transfer  $\sim 100\text{ }\mu\text{L}$  hexane layer (top) to an autosampler vial with insert for GC/MS analysis.

#### **Reference:**

1. Christie, W.W. and Han, X. *Lipid Analysis - Isolation, Separation, Identification and Lipidomic Analysis* (4th edition), 446 pages (Oily Press, Bridgwater, U.K.) (2010).