

Standard Primers

The following are available for high-throughput sequencing at no additional charge. We urge you to verify that your vector has the complimentary primer site. When submitting your samples, identify the primer you wish to use by the appropriate primer name:

Primer Name	Sequence
M13 Forward	5'-TGTA AACGACGGCCAGT-3'
M13 Reverse	5'-CAGGAAACAGCTATGACC-3'
KS	5'-CCTCGAGGTCGACGGTATCG-3'
SK	5'-CGGCCGCTCTAGA ACTAGTGGATC-3'
SP6	5'-ATTTAGGTGACACTATAG-3'
T7	5'-TAATACGACTCACTATAGGG-3'
T3	5'-ATTAACCCTCACTAAAGGGA-3'
N	5'-ACTTTATTGTCATAGTTTAGATCTATTTTG-3'
S	5'-ATAATCCTTAAAAACTCCATTTTCCACCCCT-3'
T24N (Combination of)	5'-TTTTTTTTTTTTTTTTTTTTTTTGG-3' 5'-TTTTTTTTTTTTTTTTTTTTTTTTC-3' 5'-TTTTTTTTTTTTTTTTTTTTTTTAA-3'